

EXPERIMENT 7

High-Resolution Gamma-Ray Spectroscopy

EQUIPMENT NEEDED FROM EG&G ORTEC

Source Kit SK-1G
 Other sources: 5 μCi of ^{137}Cs ; 10 μCi of ^{60}Co ; 10 μCi of ^{228}Th
 Bin and Power Supply
 GEM-10195 Coaxial Detector System (includes detector, cryostat, dewar, and preamplifier); typical specifications, 10% efficiency, 1.95-keV resolution at 1.33 MeV, 37:1 peak-to-Compton
 459 5 kV Detector Bias Supply
 572 Spectroscopy Amplifier
 444 Gated Biased Amplifier

480 Pulser
 ACE-2K MCA System including suitable IBM PC (other EG&G ORTEC MCAs may be used)
 Oscilloscope
 (Optional for Experiment 7.4) a 1- to 3-Ci americium-beryllium isotopic neutron source; if the source is not in a paraffin howitzer, place a 6-in. thickness of paraffin between the source and the HPGe detector to thermalize the source neutrons.
 ORC-7 Cable Set

Purpose

Gamma-ray energies will be measured with a high-purity germanium (HPGe) detector and research-grade electronics. The theory of response characteristics is explained, and the high-resolution results of measurement are contrasted with an NaI(Tl) scintillation detector in Experiment 3.

Introduction

Most of the experiments in this manual are written for use with EG&G ORTEC educational modules. In this experiment, which illustrates the high-resolution capabilities of HPGe detector systems, research-grade modules have been listed in order to process the pulses from the HPGe detector with a greater degree of precision to complement the detector capabilities.

Many colleges, both large and small, have Nuclear Spectroscopy Centers. In these laboratories the research efforts of the department will normally be directed in the area of high-resolution gamma-ray spectroscopy. It is possible to do a great amount of publishable research on such work as decay scheme analysis, etc., with these high-resolution systems. In many cases additional lines can be found in spectra or the improved resolution can reveal a doublet, whereas earlier measurements with NaI(Tl) detectors indicated a single energy line.

The latest decay schemes for isotopes are included in refs. 10 and 12. More recent information on certain nuclei can be obtained by writing

The Nuclear Data Project
 Building 4500
 Oak Ridge National Laboratory
 Post Office Box X
 Oak Ridge, Tennessee 37831 U.S.A.

In Experiment 3 gamma spectroscopy was studied with NaI(Tl) detectors. Typical energy resolution that can be

obtained with NaI(Tl) is $\sim 7\%$ for the 0.661-MeV ^{137}Cs gamma line. For NaI(Tl) detectors the resolution is a strong function of energy. Variation of resolution results primarily from the statistical fluctuation of the number of photoelectrons that are produced at the photocathode surface in the photomultiplier tube. Table 7.1 illustrates some typical resolutions for NaI(Tl) spectroscopy as functions of the gamma energies.

The use of germanium detectors has completely revolutionized gamma spectroscopy. Figure 7.1 illustrates the striking contrast of results obtained with these two types of

Table 7.1. Typical Resolutions of NaI(Tl) for Different Gamma Energies.

Isotope	Gamma Energy (keV)	Resolution (%)
^{166}Ho	81	16.19
^{177}Lu	113	13.5
^{133}Te	159	11.5
^{177}Lu	208	10.9
^{203}Hg	279	10.14
^{51}Cr	320	9.89
^{198}Au	411	9.21
^7Be	478	8.62
^{137}Cs	661	7.7
^{54}Mn	835	7.26
^{207}Bi	1067	6.56
^{65}Zn	1114	6.29
^{22}Na	1277	6.07
^{88}Y	1850	5.45

The information for this table was taken from *IRE Trans. Nucl. Sci.* **NS-3**(4), 57 (Nov. 1956). "Intrinsic Scintillator Resolution," by G. G. Kelley *et al.*, quoting results from F. K. McGowan, *et al.*




Fig. 7.1. ⁶⁰Co Spectrum Showing Resolutions and Peak-to-Compton Ratios for an HPGe Coaxial Detector and an NaI(Tl) Detector.

detectors. There is a factor of 30 improvement of the resolution of the data at the full-width at half-maximum (FWHM) count levels. As a result of the improved resolution, many nuclear energy levels that could not even be seen with NaI(Tl) detectors are identified easily with HPGe detectors.

In a parallel fashion the development of Si(Li) detectors has revolutionized x-ray spectroscopy. These Si(Li) devices will be studied in Experiment 8.

The purpose of this experiment is to study some of the properties of the HPGe detector systems. This experiment deals only with the practical aspects of making measurements with these detectors. To understand the properties of these detection systems, the following brief review of gamma-ray interactions and pair-production processes is included.

In Experiment 3 it was pointed out that the pair-production process at gamma energies >3 MeV is a very important gamma interaction. Figure 7.2 shows graphs for the three important gamma interactions for both germanium and silicon. The information for germanium is of interest in this experiment, and that for silicon will be used in Experiment 8.

When a gamma enters a detector, it must produce a recoil electron by one of three processes before it is recorded as an event: the photoelectric effect, the Compton effect, or pair production.

In the photoelectric process the gamma or x ray gives all of its energy to the recoil electron. It is the recoil electron that produces the electron-hole pairs in the detector that yield the output pulse. For the photoelectric process the output pulse from the detector is proportional to the energy of the




Fig. 7.2. Relative Probability of Each of the Three Types of Interactions as a Function of Energy.

gamma or x ray that produced the interaction. In the spectrum these events will show up as full-energy photopeaks.

In the Compton process there is a distribution of pulse amplitudes up to some maximum pulse height. This maximum pulse height produces the Compton edge, as explained in Experiment 3, and there is a statistical probability that each event can produce a pulse with any height up to this maximum with about an equal chance. Thus Compton events will provide a well-distributed low-energy area in the spectrum.

In modern, large detectors with high peak-to-Compton ratios, some Compton events also contribute to the full energy peak when the scattered photons undergo one or more additional interactions. This results in complete absorption.

The pair-production process can also provide a total absorption of the gamma-ray energy. The gamma enters the detector and creates an electron-positron pair. From the law of conservation of mass and energy it follows that the initial gamma must have an energy of at least 1.02 MeV because it takes that much energy to create both the negative and positive electrons. The net mass that is produced is two electron masses, and this satisfies the law of conservation of energy, $E = mc^2$.

Note that 1.02 MeV is about twice the annihilation energy that was measured from ^{22}Na (0.511 MeV). Figure 7.3 illustrates what happens in the detector in the pair-production process.

In Fig. 7.3 the e^- (ordinary electron) will produce a pulse whose magnitude is proportional to the energy of e^- (E_{e^-}). The positron will produce a pulse proportional to E_{e^+} . Since these two pulses are produced simultaneously, the output pulse from the detector would be the sum of the two pulses. When the positron annihilates in the detector, the annihilation radiation, γ_1 and γ_2 , will be produced. In Fig. 7.3 both γ_1 and γ_2 are shown escaping from the boundaries of the detector without making any further interactions. (Note: $\gamma_1 = \gamma_2 = 0.511$ MeV.) Thus, for this example, an energy of exactly 1.02 MeV escapes from the detector and is subtracted from the total energy that entered the detector. It is possible for only one, either γ_1 or γ_2 , to make a photoelectric interaction in the detector while the other escapes. In such cases the total energy absorbed is 0.511 MeV less than the original incident gamma energy. It is also possible for both




Fig. 7.3. Process of Pair Production in Germanium.

gammas to make photoelectric interactions without escaping, with all the original energy then being left in the detector. Therefore in the spectrum that is measured there will be three peaks for each gamma energy. These peaks are labeled Full-Energy Peak, Single-Escape Peak, and Double-Escape Peak, and they will be separated by 0.511-MeV increments. Figure 7.4 shows a typical spectrum that would be obtained for an incident gamma energy of 2.511 MeV. The lower end of the spectrum that includes the Compton distribution has not been included; this effect is obtained by using a biased amplifier to eliminate the lower energies and to expand the distribution of the higher-energy pulses across the range that is measured. The Single-Escape Peak occurs at $(E_\gamma - 0.511$ MeV) or 2.00 MeV, and the Double-Escape Peak occurs at $(E_\gamma - 1.02$ MeV) or 1.49 MeV. Of course, the full-energy peak represents those events for which there was a combination of pair production and photoelectric effect in which all the energy was absorbed in the detector.




Fig. 7.4. Typical Spectrum for an Incident Gamma Energy of 2.511 MeV.

Now refer again to Fig. 7.2 and specifically to the curves for the interactions in germanium. The absorption cross section, plotted in the y direction, is a measure of the relative probability that an interaction will take place in a germanium detector. These probabilities of relative interactions, for the most part, determine the shape of the observed spectrum. For example, a photon (or gamma) with an energy of 100 keV has an absorption cross section of ~ 55 barns/atom for the photoelectric process. The corresponding Compton cross section is ~ 18 barns/atom. There is no pair production. This indicates that at 100 keV there are 3 times as many photoelectric interactions as Compton interactions, since this is the approximate ratio of the cross section. Figure 7.5 shows the shape of a spectrum that could be expected for measurement of the 100-keV energy events.

The sum of counts under the photopeak would be 3 times the sum under the Compton distribution. For larger crystals Σ_{pp} would be even >3 times Σ_c because some of the scattered gammas from the Compton interactions would make photoelectric interactions before escaping from the crystal. For an infinitely large crystal there would be no Compton distribution since the crystal would then totally absorb all of the incident gammas.




Fig. 7.5. Typical Spectrum Expected for 100-keV Energy in HPGe.

The shape of the spectrum changes drastically from 100 keV to 1 MeV. Figure 7.6 shows the gamma spectrum that could be expected for the 1-MeV gammas incident on an HPGe detector. From Fig. 7.2 the ratio of Compton cross section to photoelectric cross section is ~ 90 ; so in Fig. 7.6, Σ_c is 90,000 and Σ_{pp} is 1000. The variation of cross sections for HPGe and Si(Li) detectors can also be approximated from Fig. 7.2. For example, at $E_\gamma = 400$ keV the photoelectric cross section for germanium is 6 barns/atom and that for silicon is ~ 0.1 barn/atom. This is a ratio of 60:1 and indicates that there will be 60 times as many counts under the photopeak for a germanium detector as for a silicon detector at 400 keV, assuming that the detectors are the same size. The reason for this is that the photoelectric cross section varies as Z^5 , where Z is the atomic number of the absorbing material. The atomic number of Ge is 32 and is 14 for Si. The ratio of these two numbers raised to the 5th power is 62.22, and this is within remarkable agreement with the above cross-section ratios.




Fig. 7.6. Typical Spectrum Expected for 1-MeV Energy in HPGe.

EXPERIMENT 7.1

Energy Resolution with an HPGe Detector

The instructor will provide the HPGe detector and the instructions for its use. Read the instruction manual carefully before attempting to use the detector. This is a very expensive detector system and must be handled carefully.

Procedure

1. Install the 459, 480, 572, and 444 in the bin and power supply and interconnect the modules as shown in Fig. 7.7. The preamplifier will be mounted on the HPGe detector and the interconnection for the signal connections to the detector is made through the preamplifier. Connect the 459 to the Detector Bias Input and connect the 480 to the Test Input on the preamplifier. Connect the power cable for the preamplifier to the Preamplifier Power connector on the 572. Connect a signal output from the preamplifier to the 572 input. Connect the Unipolar Out from the 572 to the linear input of the 444, and connect the 444 output to the input of the MCA. Set the module controls as follows:

572 Amplifier: Positive input (verify with instructor); Unipolar Output; Shaping time 6 μ s; Delay Out.

444 Gated Biased Amplifier: Coarse Gain X2; Fine Gain 1/2; Bias Level 20/1000; Normal mode; DC-Couple; Pulse Duration 6 μ s (internal control); Anticoincidence and Internal Strobe (rear-panel controls).

480 Pulser: Attenuated output.

459 5 kV Detector Bias Supply: Leave at zero until all other connections and adjustments have been made; consult the instructions for the HPGe detector to determine both the polarity and the amplitude of bias that are to be used, and apply the correct amount of bias in the correct polarity when ready to operate.

2. Place the ^{60}Co source from source kit SK-1G ~ 1 cm from the face of the detector. Adjust the gain of the 572 Amplifier so that the 1.33-MeV gamma has an amplitude of 6 V at the amplifier output. The two lines for 1.17 and 1.33 MeV should be quite easily seen on the oscilloscope. Check the output from the 444 Gated Biased Amplifier to make sure it looks reasonable.

3. Observe the display of the spectrum on the MCA. Adjust the bias level and gain on the 444 and the gain on the 572 until the two sharp photopeaks are positioned as shown in Fig. 7.8. In this measurement the two photopeaks should be separated by at least 200 channels, based on 1024 channels total.

4. From the positions of the two photopeaks make a calibration curve of energy (y direction) vs channel number (x direction) and determine the keV/channel.

EXERCISES

- a. What is the resolution in keV for the two photopeaks? How does this compare in value with the detector's resolution specifications?

- b. From the data, determine the energies of the Compton edges for the two gammas. How do these compare with the values that were calculated from the formula used in Experiment 3?




Fig. 7.7. Electronic Interconnections for Experiment 7.1.

5. Turn on the 480 Pulser and adjust the output so that the pulser peak falls about half-way between the 1.17- and 1.33-MeV peaks. After the pulser has produced enough counts for its peak channel to have 1000 counts, read that portion of the analyzer memory and determine the resolution of the pulser. This is the electronic resolution, $R(E)$. The contribution from the detector to the overall resolution can be calculated from the formula

$$\text{system resolution} = \sqrt{[R(d)]^2 + [R(E)]^2} \quad (1)$$

where $R(d)$ is the detector resolution and $R(E)$ is the electronic resolution. These resolutions are said to add in quadrature. There is a lower limit to $R(d)$ which is energy-dependent. The recoil electron that is produced in the gamma interaction loses energy in the HPGe detector by dE/dx . The average energy required to produce an ionization in germanium is 2.95 eV/electron-hole pair. Thus for a 1.5-MeV recoil electron there would be 5.08×10^5 electron-hole pairs produced. The production of electron-hole pairs is a process that is statistical in nature, and hence there are fluctuations in the actual number produced. When the proper statistics are used, the theoretical lower limit to $R(d)$ is given by

$$R(d) = K \sqrt{F \cdot E} \quad (2)$$




Fig. 7.8. Locating the Two Photopeaks for ^{60}Co in the HPGe Spectrum.

where K is a constant, E is the energy of the photon in MeV, and F is the statistical Fano factor. To a very good approximation this equation reduces to

$$R(d) \text{ (in keV)} = 1.35 \sqrt{E \text{ (in MeV)}}. \quad (3)$$

Solving Eq. (3), the theoretical lower limit of detector resolution is 1.44 keV for a 1-MeV gamma and is 4.5 keV for a 10-MeV gamma.

EXERCISES

c. Calculate the values of $R(d)$ from Eq. (3) for the values of E in Table 7.2.

Table 7.2

Energy (MeV)	Theoretical $R(d)$ (keV)
0.1	
0.3	
0.5	
1.0	
3.0	
6.0	
8.0	
10.0	

d. Make a plot of the values from Table 7.2 on linear graph paper. From Eq. (1) calculate the experimental $R(d)$ for the 1.33-MeV peak of ^{60}Co . How does this compare with the theoretical limit? Remember that the $R(d)$ theory is the absolute lower limit of the resolution value.

EXPERIMENT 7.2

Photopeak Efficiency for HPGe Detectors

Resolution with HPGe detectors is better by a factor of 30 or more than that obtained with NaI(Tl) conventional detectors. Coupled with this dramatic increase in resolution is a compromise of the photopeak efficiency. The pricing of HPGe detectors is related to their photopeak efficiency. The standard method for comparing the efficiencies of HPGe detectors with NaI(Tl) detectors is to compare their counting rates at the 1.33-MeV line of ⁶⁰Co, using a standard distance of 25 cm from the source to the detector face and placing the source on the detector axis.

The resolution of the HPGe detectors is so many times better than that of the NaI(Tl) that the ability to see a photopeak above the Compton distribution is considerably enhanced. Consider a simple example in which the efficiencies of the HPGe and NaI(Tl) are assumed to be the same. In a particular experiment we observe 10,000 counts under the photopeak for each detector. If the resolution of the HPGe detector is only 10 times that of the NaI(Tl) detector, the HPGe photopeak will have 10 times the maximum number of counts that the NaI(Tl) detector has, because the area under the photopeak (10,000 for this example) is approximately proportional to the width times the height of the peak. Since the width of the HPGe peak is only 1/10 the width of the NaI(Tl) peak, its height must be 10 times as great.

This example can easily be extrapolated to real situations where the advantages of superior resolution are very important. For example, Fig. 7.9 shows the striking differences for a spectrum obtained on a mixed sample of ⁷⁶As, ¹²²Sb, and ¹²⁴Sb with each of the two types of detectors. Each of the closely spaced energy lines is shown separately in the HPGe spectrum and are all included in a single broad photopeak in the NaI(Tl) spectrum.

In this experiment, we will measure some of these photopeak efficiencies and also determine the peak-to-Compton ratio for an HPGe detector.

Procedure

1. Use the same equipment setup that was used for Experiment 7.1. Adjust the gain of the 572 and bias level and gain of the 444 for a ⁶⁰Co spectrum similar to that of Fig. 7.10.
2. Accumulate the spectrum in the MCA for a time period long enough to determine heights h_1 and h_2 to a fair degree of accuracy. In Fig. 7.10, h_1 is the 1.33-MeV photopeak and h_2 is the maximum for the comparable Compton distribution, normally located just below the Compton edge. Read the data from the analyzer.




Fig. 7.9. Comparative Spectra Taken with HPGe and NaI(Tl) Detectors.




Fig. 7.10. Typical Distribution of ⁶⁰Co Spectrum for Experiment 7.2.

EXERCISE

a. Calculate the peak-to-Compton ratio, which is h_1 divided by h_2 in Fig. 7.10. Compare your value with the value for this ratio for the detector; check with your laboratory instructor for the record of the ratio.

3. Clear the spectrum from the MCA. Place a $10\text{-}\mu\text{Ci } ^{60}\text{Co}$ source at a distance of exactly 25.0 cm from the face of the detector.

4. Accumulate a spectrum for this source for a time period long enough for the sum of counts under the 1.33-MeV photopeak to be about 3000. Read the data from the analyzer and be sure to record the live time for the measurement.

EXERCISES

b. Calculate the number of counts per second for the events that were recorded in the 1.33-MeV photopeak; call this R_1 :

$$R_1 = \frac{\Sigma_{pp}}{\text{time in seconds}} \quad (4)$$

c. The rate, R_1 , from Eq. (4) is to be compared with the rate, R_2 , that is expected for the same source when it is located 25.0 cm from the face of a 3- x 3-in. NaI(Tl) detector. The efficiency of this size NaI(Tl) detector for a source-to-detector distance of 25.0 cm is given as 1.2×10^{-3} , which is from the "Gamma Ray Spectrum Catalog," by R. L. Heath, Idaho Falls Report IDO-16880. Using ϵ_1 for this number, the number of counts, (N), that you would observe under the photopeak for a 3- x 3-in. NaI(Tl) detector at 25.0 cm source distance is given by

$$N = \epsilon_1 A t, \quad (5)$$

where A is the gamma activity of the source in counts per second and t is the live time in seconds. The rate R_2 will then be

$$R_2 = \frac{N}{t} = \epsilon_1 A. \quad (6)$$

Since ^{60}Co has a 1.33-MeV gamma ray for each decay, A is given by

$$A = 3.7 \times 10^4 (x), \quad (7)$$

where x is the source strength in microcuries (μCi).

Calculate R_2 from Eq. (6). The relative photopeak efficiency is obtained for the detector by

$$\text{relative photopeak efficiency} = \frac{R_1}{R_2} \times 100. \quad (8)$$

Calculate this value for your measurement and compare it with the value that is recorded for the detector; check with your laboratory instructor for the record of the detector's efficiency.

EXPERIMENT 7.3

Escape Peaks and Efficiency for HPGe Detectors

As discussed earlier, when an incident gamma with sufficient energy enters the crystal it can create an electron-positron pair. When the positron annihilates, two gammas with equal energy at 0.511 MeV are produced and these leave with an angular separation of 180° . In Fig. 7.3 these two gammas are shown as γ_1 and γ_2 . For small detectors it is very probable that both γ_1 and γ_2 will escape from the detector before they make any further interactions in the crystal. The energy thus absorbed would be $E_\gamma - 1.02$ MeV and is shown as the Double-Escape Peak in Fig. 7.4. As the detector size is increased, the probability is greater that either γ_1 or γ_2 will make a photoelectric interaction within the crystal. If one of these gammas does make a photoelectric interaction, the energy of the event that is recorded in the detector is the Single-Escape Peak in Fig. 7.4. For even larger detectors the probability of photoelectric interactions is further increased when both γ_1 and γ_2 interact and the total energy of the gamma is absorbed in the crystal. Figure 7.11 shows some measurements that have been made for coaxial and planar HPGe detectors. From this figure the ratios of Full-Energy, Double-Escape Peak, and Single-Escape Peak efficiencies can be determined by inspection for the size of detector that is identified in the figure.

To see how the measurements were made for Fig. 7.11, consider the E_γ of 2.511 MeV shown in Fig. 7.4. Assume that the




Fig. 7.11. Measured Efficiencies for HPGe Detector.




Fig. 7.12. Typical High-Energy Gamma Spectrum from a Neutron Source with Iron Scatterer.




Fig. 7.13. Typical Spectral Response of an EG&G ORTEC HPGe Detector to Various Gamma-Ray Energies.

Compton distribution has been subtracted for each peak in Fig. 7.4 and that the following sums have been measured:

- Σ at Full Energy, 2.511 MeV = 6000,
- Σ at Single-Escape, 2.00 MeV = 1000,
- Σ at Double-Escape, 1.489 MeV = 3000.

From these numbers the simple ratios can be obtained.

Procedure

1. Use the same equipment setup that was used for Experiment 7.1. Use the ^{60}Co and ^{137}Cs gamma sources from SK-1G. Adjust the system gain and bias on the 572 and 444 to calibrate the analyzer roughly from 1 to 3 MeV.

2. Remove the energy calibration sources and use a ^{228}Th (or other high-energy) source to accumulate a spectrum. Accumulate for a period of time long enough to see all the pronounced peaks in the spectrum. Read the data from the analyzer.

EXERCISES

a. Plot the spectrum on semilog graph paper. On the plot identify all the major peaks and the corresponding escape peaks. Compare the energies at these peaks with those that are identified with the source.


b. Calculate the escape-peak ratios. Define Σ_f as the sum of counts under the full-energy peak, Σ_1 as the sum under the single-escape peak, and Σ_2 as the sum under the double-escape peak. Be sure to subtract the Compton distribution from these sums. Then determine the ratios Σ_f/Σ_1 , Σ_f/Σ_2 , and Σ_2/Σ_1 . How do these ratios compare with those the laboratory instructor has for the ^{228}Th source and the detector you are using?

EXPERIMENT 7.4

(Optional, recommended if Experiments 16, 17, or 18 are to be done)

The Response of HPGe Detectors to High-Energy Gammas

If an isotopic neutron source of the Am-Be type is available, it is possible to obtain high-energy gammas from this source. The neutrons from the source are produced by



The Q value for the reaction is ~ 5 MeV. Since the alpha energies from most sources are also ~ 5 MeV, it is possible to produce neutrons with these sources up to 10 MeV. The neutron spectrum from these sources shows a distribution of neutron energies up to this maximum energy of ~ 10 MeV. What is of more importance is that in the reaction it is also possible for ^{12}C to be left in an excited state. The de-excitation of ^{12}C is by gamma emission. Gammas from the second excited state of ^{12}C have an energy of 7.656 MeV and make an excellent source of high-energy gammas.

Procedure

1. Use the same equipment setup that was used for Experiment 7.1. Use the ^{60}Co source and the pulse generator to adjust the gain and bias of the system so that the MCA range is ~ 3 to 8 MeV.

2. Use a block of paraffin ~ 6 in. thick between the detector and the neutron source, and place the source ~ 12 in. from the detector. The paraffin will thermalize the neutrons from the source without attenuating the high-energy gammas. In some cases the neutron source is in a paraffin howitzer; if this is the case, place the source close to the outside of the howitzer.

3. Accumulate a spectrum. This will require several hours, and sometimes overnight runs are necessary. Read the data from the analyzer.

EXERCISE

Plot the spectrum on semilog paper and identify all the peaks. As in Experiment 7.3, calculate the ratios Σ_f/Σ_1 , Σ_f/Σ_2 , and Σ_2/Σ_1 .

For your reference, Fig. 7.12 is a plot on linear graph paper for a typical neutron source with an iron scatterer. The reaction is ${}^{56}\text{Fe}(n,\gamma){}^{57}\text{Fe}$, which yields two high-energy gammas from the ${}^{56}\text{Fe}$ scatterer. The main photopeak energies are 7.639 and 7.625 MeV.

The resolution of a gamma line is dependent on the gamma energy of the peak. The student can easily verify this by measuring the resolution of the detector for the sources available in the laboratory. Figure 7.13 shows the spectral response and resolution of several common sources for an HPGe detector.

References

1. A. Coche and P. Siffert, *Lithium Drifted Silicon and Germanium Detectors in Semiconductor Detectors*, G. Bertolini and A. Coche Eds., Elsevier-North Holland, Amsterdam (1968).
2. J. C. Phillpott, *IEEE Trans. Nucl. Sci.* **NS-17**(3), 446 (1970).
3. F. Adams and R. Dams, *Applied Gamma Ray Spectrometry*, 2nd Edition, Pergamon Press, Ltd. (1970).
4. U. E. P. Berg, H. Wolf, et al., *Nucl. Instrum. Methods* **129**, 155 (1975).
5. C. Meixner, *Nucl. Instrum. Methods* **119**, 521 (1974).
6. D. M. Walker and J. M. Palms, *IEEE Trans. Nucl. Sci.* **NS-17**(3), 296 (1970).
7. M. L. Stelts and J. C. Browne, *Nucl. Instrum. Methods* **133**, 35 (1976).
8. 20th Nuclear Science Symposium, *IEEE Trans. Nucl. Sci.* **NS-21**(1) (1974).
9. 19th Nuclear Science Symposium, *IEEE Trans. Nucl. Sci.* **NS-20**(1) (1973).
10. C. M. Lederer and V. S. Shirley, Eds., *Table of Isotopes*, 7th Edition, John Wiley and Sons, Inc., New York (1978).
11. G. F. Knoll, *Radiation Detection and Measurement*, Wiley, New York (1979).
12. Nuclear Data Sheets, Edited by the Nuclear Data Project, Academic Press, New York and London.